



Release Notes

QuBAS 3.1.5

FileHold ID: 3318

FH Version: 1

Approved: See metadata

Page 1 of 2



Version / Release	3.1.5 / 08 October 2025
3.1.5 Introduction	3.1.5 is a minor release that fixes a couple of issues with web browsers.
3.1.5 Upgrade advice and versions supported	3.1.4 is a recommended but not required upgrade from 3.1.y. It contains a couple of fixes for browser issues. The analysis engine is identical to that of version 3.1.4. Previous results are not impacted. From this release, the supported versions of QuBAS are: 3.0.1, 3.1.0, 3.1.2, 3.1.4 If you require help or advice about upgrade please contact support@quantics.co.uk
3.1.5 New features	This version is dedicated to bug fixes and doesn't have new features.
3.1.5 Bug fixes	❑ Chrome users don't get prematurely logged out. ❑ Leaving QuBAS open for a long time (~24 hours) doesn't take up a lot of memory. ❑ Some combination of special characters in the password of the email server were not allowed. This is no longer the case.
3.1.5 Known Issues	❑ The logging out message from the development room only disappears with a screen refresh. ❑ QuBAS fails to load data files with a path name >260 characters due to Windows default settings. Keep filenames and path to < 260 characters.
3.1.5 Testing	The analysis engine has not changed since version 3.1.4 and thus didn't require further testing.
3.1.5 Build	Build #3
3.1.5 Installer	The installer is QuBAS_3.1.5_3.exe The MD5 is: ED4C9240625577CD2852570BAE8F9F34 QuBAS_3.1.5_3.exe